



MICROARRAY TERMINOLOG

- Feature--an array element
- Probe--a feature corresponding to a defined sequence
- Target--a pool of nucleic acids of unknown sequence

Cancer Genetics Branch





Kinds of array elements

- Synthetic Oligonucleotides
- PCR products from Cloned DNAs

Genomic DNA

Cloned DNA







Microarray Manufacture

- Printing
- Synthesis in situ


Cancer Genetics Branch





MICROARRAY


Retended Tuantity of target bound to each probe in a complex hybridization

- Must have high sensitivity, low background
- High spatial resolution essential
- Dual channel capability preferred











SINGLE NUCLEOTIDE POLYMORPHISMS

- •Polymorphic SNPs occur approximately every 1 kb in the human genome.
- •Dense SNP maps provide a basis to design microarrays for genome scanning

Cancer Genetics Branch



LABELLING SNPs

Genomic DNA




Unlabeled amplicons




Labeled amplicons



Hybridize to microarray





National Human Genome Research Institute


ACCURACY OF SNP CHIP


Table 3. ABACUS SNP Detection and Genotyping Accuracy

A. Accuracy of autosomal SNPs detection	Mandelland	Total Decilate
Singleton SNPs	Verified 1.7	Total Possible
Non-singleton SNPs	91	91
Total SNPs	108	108
	100	100
B. Number of autosomal SNPs electronically verified		
Number of SNPs electronically verified	371	
C. Accuracy of autosomal genotype calls		
Number of verified homozygous genotype calls	1515	
Number of incorrect homozygous genotype calls	0	
Percent correct homozygote calls	100.00%	
Number of verified heterozygous genotype calls	423	
Number of incorrect heterozygous genotype calls	3	
Percent correct heterozygoté calls	99.30%	
D. Accuracy of haploid genotype calls		
Number of bases sequenced (6X coverage)	17,423	
Number of bases different from microarray chip calls	0	
Percent of bases identical	100.00%	

Cutler DJ et al. Genome Res. 2001 11:1913-25.

Cancer Genetics Branch





DNA Microarray Applications

Resequencing

Mutations Polymorphisms


Gene copy number


Cancer Genetics Branch





Gene Copy Number


- Array format CGH
- Large insert clones
- cDNA clones/exons
















DNA Microarray Applications

Resequencing

Mutations Polymorphisms


- Gene copy number
- Gene expression


Cancer Genetics Branch



High throughput analysis of gene expression

- cDNA library sequencing
- Serial analysis of gene expression (SAGE
- Microarray hybridization







Oligo versus cDNA Arrays for Expression Analysis

Cancer Genetics Branch

National Human Genome Research

Oligonucleotide Arrays: Pros

- Complete control over sequence
- Sequence and geometric perfection
- Extremely high feature density



Oligonucleotide Arrays: Cons


- Lack of Flexibility in Some Formats
- Absolute Requirement for Sequence Da
 - •Risk of uneven Performance by Individual Array Elements (Lack of Oligo Picking Rules)


Cancer Genetics Branch




cDNA Arrays: Pros


- High Degree of Flexibility
- Sequence Independent
- High Stringency Hybridization
 - High Signal Intensity: No Need for Signal Amplification




























Clustering of Human ESTs

Final Number of Clusters (sets) (Unigene 146)

96109 sets total

21857 sets contain at least one known gene

94916 sets contain at least one EST

20664 sets contain both genes and ESTs

193 genes are not represented by an identifiable EST.

Cancer Genetics Branch



Histogram of cluster sizes for UniGene build 146

	Cluster	Number
	size	of clusters
	1	33700
	2	13467
	3-4	15267
	5-8	10197
06 000 alvetora	9-16	5777
96,000 clusters	17-32	3894
	33-64	3549
	65-128	4031
	129-256	3718
	257-512	1732
	513-1024	537
	1025-2048	162
	2049-4096	56
	4097-8192	19
	8193-1638	34 3

Cancer Genetics Branch


Expressed Sequence Tags:
Options for Array Constructi

*"Standard" clone sets

Custom clone sets

Synthetic oligonucleotides

Cancer Genetics Branch





APPLICATIONS OF EXPRESSION ARRAYS

Direct comparisons (Induction)

Biological system critical

Expression profiling

Requires statistical tools

Power arises from increasing sample number


Cancer Genetics Branch



APPLICATIONS OF EXPRESSION ARRAYS

Statistical tools for large datasets


First generation approaches.



APPLICATIONS OF EXPRESSION ARRAYS:

- TБИМОВ PROFILING
 taxonomy of cancer
- Methods lead to gene identification
- Individualized diagnosis and therapy

Cancer Genetics Branch



APPLICATIONS OF EXPRESSIONARRAYS: GENE IDENTIFICATION

- Groups of genes
 - Pathways
 - Co-regulated
 - Correlate with copy #
 - Correlate clinically
- Candidate disease genes

Cancer Genetics Branch



APPLICATIONS OF EXPRESSION ARRAYS: TUMOBLER FILING

Unsupervised
Supervised
•Classification

 Can classify with respect to any clinically interesting variable

Cancer Genetics Branch

Alveolar Rhabdomyosarcoma


Pax3 chromosome 2


Method

- •Compared 7 ARMS with 6 unrelated cancers cell lines
- •Using cDNA microarray containing 1238 elements

Cell Line Characteristics


Cell Line	Pax3-FKHR	Diagnosis
ARMS1	1	ARMS
RH3	+	ARMS
RH4	+	ARMS
RH5	+	ARMS
RMS18	+	ARMS
RMS13	+	ARMS
RH28	+	ARMS
A204	-	Undifferentiated Sarcoma
NGP127	-	Neuroblastoma
TC71	-	Ewing's Sarcoma
UACC-903	-	Melanoma
PC3	-	Prostate Carcinoma
MDA-MB-4	36 -	Breast Carcinoma
Control		
NIL-C	-	Fibroblast







Matrix of Pearson Correlation Coefficients Distance Map


78 pair-wise comparisons





Hierarchical Clustering Dendrogram























Molecular Taxonomy of Small Round Blue Cell Tumors Hypothesis

 Using cDNA microarrays we can identify the genes whose expression level is

• Utilize these genes to classify the characteristic for that cancer & type small blue round cell tumors into the

correct disapportis sotogories



cDNA Microarray Analysis

Experiments

Burkitt's Lymphoma 8

EWS-Tumor

13

EWS-Cell line

10

Neuroblastoma

12

>500, 000 Data Points

RMS-Tumor

10

RMSACHILIDE TRAY

4,000 sequence verified known genes


2,55t/Unknownce verified EST


2,567 genes


Artificial Neural Networks


Pattern Recognition


Training


























WHAT WE HOPE TO LEARN IN THE FUTURE

- IMPROVE THE DIAGNOSTIC CATEGORIZATION OF TUMORS.
- IDENTIFY USEFUL PREDICTIVE MARKERS FOR OUTCOME AND THERAPEUTIC RESPONSE (ARRAY OR CONVENTIONAL).
- IDENTIFY POINTS FOR INTERVENTION:

CRITICAL PATHWAYS

DRUG TARGETS

Cancer Genetics Branch



CLINICAL CORRELATIVE STUDIES USING MICROARRAYS

- DEFINE QUESTION AND PATIENT SAMPLE.
- APPROPRIATE AND RIGOROUS STATISTICAL ANALYSIS OF ARRAY DATA.
 - RESULT: GENES WHICH CARRY INFORMATION RELEVANT TO QUESTION POSED.
 - DEVELOP FORMAL CLASSIFIER.
- VALIDATE ON ADDITIONAL SAMPLE SET.



MODEL SYSTEM WITH CLEAR THERAPEUTIC IMPLICATIONS: GASTROINTESTINAL STROMAL TUMOR

- RELATED TO THE INTERSTITIAL CELLS OF CAJAL
- KIT MUTATIONS
- STI-571 SENSITIVITY
- THE BEST "CREDENTIALED" TARGETS ARE THOSE ACTIVATED BY MUTATION.


Cancer Genetics Branch




SAMPLES

- 13 MALIGNANT GISTs
- ALL KIT POSITIVE BY IHC
- ALL WITH KIT MUTATIONS
- 4 GUT WALL PRIMARIES
- 8 INTRA-ABDOMINAL EXTENSION
- 1 LIVER METASTASIS
- 6 COMPARISON TUMORS: EXTRA-GI SPINDLE CELL MORPHOLOGY

Cancer Genetics Branch









DATA ANALYSIS

- PREFILTER FOR QUALITY AND IDENTIFY GENES HIGHLY EXPRESSED IN GIST: 1987 GENES
- RANK BY WEIGHTED DISCRIMINATOR METHOD $w(g, \pm) = \mu_+(g) \mu_-(g) / [_+(g) + _-(g)]$
- RANDOM PERMUTATION TEST (105 trials)
- CLUSTER ANALYSIS USING DISCRIMINATORS

Cancer Genetics Branch



Mational Human		TOP DIS	CRIMINA	ATORS FOR GIST
Research	Rank	<u>Weight</u>	Alpha	Gene Description
and the second		7.55575	O	v-kit sarcoma oncogene
10/1/25	2	6.48306	0	G coupled receptor 20
	3	4.60057	O	G coupled receptor 20
0,00	4	4.51681	O	annexin A3
	5	3.33057	O	KIAA0353 protein
	6	3.31734	O	phosphofructokinase, muscle
	7	2.95095	0.000	008 DKFZP434N161 protein
	8	2.83435	O	protein kinase C, theta
	9	2.79721	O	butyrylcholinesterase
1	10	2.72752	O	annexin A3
				Cancer Genetics Branch






CONCLUSIONS

- MALIGNANT GISTS EXHIBIT A DISTINCT AND HIGHLY COHERENT GENE EXPRESSION PROFILE.
- THIS GENE LIST IS RELEVANT BOTH TO GIST GROWTH AND NORMAL ICC FUNCTION.
- KIT, A CRITICAL GENE IN REGULATING GIST GROWTH, IS THE BEST DISCRIMINATOR FOR THIS DISEASE.
- EXTENDING THIS APPROACH TO OTHER CANCERS MAY HELP IDENTIFY NEW DISEASE SPECIFIC DRUG TARGETS.

Cancer Genetics Branch



A RECURRING PROBLEM

Oncogenes

Transcription factors

Hormones/growth factors

Drugs

Toxins

Radiation

Downstream Genes

- Direct targets
- Indirect targets

Cancer Genetics Branch

Retroviru S: Empty vestor Paxs Paxs Fikhr


RESULTS

Mouse cDNA array 2200 genes

PAX3-FKHR vs Empty Vector

PAX3
vs
Empty Vector

3T3 PARENT vs
Empty Vector



- 1. Troponin C
- **2. IGFBP5**
- 3. Myogenin
- 4. Six1
- 5. Troponin T
- 6. IGF2

